



INCOR® - New Intelligent Charging Unit

In the context of our continuous product improvements we adapted the INCOR® charging unit. You and your patients can now easier handle the charging and discharging process of the batteries.

The new intelligent charging unit is part of the INCOR® system and charges as well as calibrates the batteries. This product improvement enables you to reactivate completely discharged batteries. Furthermore, with the new LED displays you have a clear overview of the batteries' charging status.



What is new?

- Distinct icons improve the overview on the charging status of the batteries
- The charging unit automatically identifies and informs the user if calibration of the batteries is required
- The user decides when the best point in time has come to start the calibration process



	LED	Icon	Color
1		Battery with lightning	green/ orange
2		Wrench	orange
3		Warning triangle with exclamation mark	red

The exact meaning is detailed in the instructions for use

Benefits of the new charging unit

- The charging process is optimized for each battery
- More safety and mobility for the patient due to a clear overview of the batteries' charging status
- The batteries have a longer operation life with regular calibrations
- Even completely discharged batteries can normally be reactivated by the user
- Generally easier handling and improved overview